
# Roadability

ENGLISH	3 - 5
DEUTSCH	6 - 8
NEDERLANDS	9 - 11
DANSK	12 - 14
NORSK	15 - 17
SVENSKA	18 - 20
SUOMI	21 - 23
ESPAÑOL	24 - 26
FRANÇAIS	27 - 29
ITALIANO	30 - 32
PORTUGUÊS	33 - 35
PORTUGUÊS DO BRASI	36 - 38
РОССИЯ	39 - 41
POLSKI	42 - 44
ČESKY	45 - 47
SLOVENSKY	48 - 50
中國	51 - 53
日本	54 - 56
ΕΛΛΗΝΙΚΗ	57 - 59

## ENGLISH


## Roadability

Driving UP kerbs (max. 100 mm)

## Forwards


This technique is recommended for the experienced users.

- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb
- Balance on the rear wheels and lift the front wheels from the ground. Bend forward and pull hard on the hand rims.

## Forwards, with assistance


- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb
- Tilt the chair backwards to lift the front wheels above the kerb.
- Lift on the handles to drive the rear wheels above the kerb.

## Backwards

This technique is only to be used at low edges and is depended on the height of the footplate.

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb
- Bend forward and pull hard on the hand rims.

## Backwards, with assistance

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb.
- Tilt the chair backwards to lift the front wheels.
- Lift on the handles to pull the chair up. Pull the chair away from the kerb and put down the front wheels.


## ENGLISH

EN

### Driving DOWN kerbs (max. 10 cm)

#### Forwards

- This technique is recommended for the experienced users.
- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb
- Pull hard on the hand rims and drive forward. All wheels will touch the ground at the same time.

#### Forwards, with assistance


- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb.
- Tilt the chair backwards to lift the front wheels.
- Drive carefully down the kerb and put down the front wheels.

#### Backwards

This technique must never be used if the height difference is more than 10 cm and only by experienced users.

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb
- Bend forward and drive carefully down the kerb.
- High tilting risk at this technique

#### Backwards, with assistance

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb.
- Drive carefully down the kerb and pull the chair away from the kerb.
- Put down the front wheels.

## ENGLISH

## Driving UP stairs

- Always ask for help
- Never enter an escalator even with help from others

## Backwards, with assistance

- Make sure the anti-tip device is deactivated and the handles is tightened securely
- Drive backwards toward the first step
- Tilt the chair and pull up the chair one step at a time. Make sure to keep balance on the rear wheels
- Pull the chair away from the last step and put down the front wheels

## Driving DOWN stairs

- Always ask for help
- Never enter an escalator even with help from others

## Forward, with assistance

- Make sure the anti-tip device is deactivated
- Drive towards the first step and tilt the chair
- Drive down one step at a time. Make sure to keep balance on the rear wheels
- Put down the front wheels after the last step
- If two helpers are present, one can assist by lifting at the front
- The helpers need to make sure to lift the correct way - with the strength from the legs and not from the back

## Driving DOWN/UP a hill

Follow the below to avoid tilting when driving down/up a hill.


## Driving down:

Sit upright in the chair. Correct the speed with the hand rims and not with the brakes.


## Driving up:

Bend forward to correct the point of gravity.

- Avoid to turn around in the middle of the hill.
- Always sit as upright as possible
- Better to ask for help than to take a risk on your own.

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## R82 Inc.
## R82 au Pty Limited